

1105 North Tennessee Street Cartersville, Georgia 30120

FAX 770-382-4757 www.seanimallabs.com

BOVINE PREGNANCY TESTING INSTRUCTIONS

Preparation:

- 1. Confirm all cows to be tested are 28 days post-breeding and 90 days post-calving.
- 2. Collect needles, tubes and holders.
- 3. Number vials sequentially with a tube number beginning with 1.

Testing:

- 1. Collect a 2 cc blood sample beginning with tube number 1.
- 2. After the tube number, write cow identification number on the tube.
- 3. Use a new needle for each cow.





Packaging and Shipping:

- 1. Complete a sample submission form. Write the animal ID on the submission corresponding to each tube number.
- 2. If you are not shipping samples the same day refrigerate the samples overnight.
- 3. Package the samples for shipping with the submission form (no ice is necessary).
- 4. Ship to: Southeastern Animal Labs 1105 North Tennessee Street

Cartersville, GA 30120